- EasyT -
- OPN NC -
- BEO NC -
- NIC Nano *hydro* -
- NIC 1.1 *hydro* -

PTCA Dilatation Catheters

Language				Page
us	English (US)	Instructions for Use	PTCA Dilatation Catheter	8
en	English (Non-US)	Instructions for Use	PTCA Dilatation Catheter	10
de	Deutsch	Gebrauchsanweisung	PTCA Dilatationskatheter	12
fr	Français	Mode d'emploi	Cathéter de dilatation à ballonnet pour ACTP	14
bg	Български	Инструкции за употреба	Дилатационен катетър за ПТКА	16
zh	简体中文	使用说明	PTCA 扩张导管	18
cs	Česky	Návod k použití	Dilatační katétr pro PTCA	20
da	Dansk	Brugsanvisning	PTCA-dilatationskateter	22
et	Eesti keel	Kasutusjuhend	PTKA dilatatsioonikateeter	24
es	Español	Instrucciones de uso	Catéter de dilatación para ACTP	26
fi	Suomi	Käyttöopas	Sepelvaltimon Pallolaajennuskatetri	28
el	Ελληνικά	Οδηγίες χρήσης	Καθετήρας διαστολής για PTCA	30
hr	Hrvatski	Upute za upotrebu	Dilatacijski kateter za PTCA	32
hu	Magyar	Használati útmutató	PTCA tágitókatéter	34
it	Italiano	Istruzioni per l'uso	Catetere a palloncino per PTCA	36
ko	한국어	사용 설명서	PTCA 확장용 카테터	38
lt	Lietuvių k.	Naudojimo instrukcijos	PTKA plėtimo kateteris	40
lv	Latviski	Lietošanas pamācība	PTCA dilatācijas katetrs	42
mt	Malti	Istruzzjonijiet għall-Użu	Kateter għal Dilatazzjoni f'PTCA	44
nl	Nederlands	Gebruiksaanwijzing	PTCA-dilatatiekatheter	46
no	Norsk	Bruksanvisning	PTCA-dilatasjonskateter	48
pl	Polski	Instrukcja obsługi	Cewnik rozszerzający do PTCA	50
pt	Português	Instruções de Utilização	Cateter de dilatação PTCA	52
ro	Română	Instrucțiuni de utilizare	Cateter de dilatare pentru ACTP	54
ru	Русский	Инструкция по применению	Дилатационный катетер для ЧТКА	56
sr	Srpski	Uputstvo za upotrebu	Balon dilatacioni kateter	58
sv	Svenska	Bruksanvisning	PCI-ballongkateter	60
sl	Slovenščina	Navodila za uporabo	Dilatacijski kateter za PTCA	62
sk	Slovenský jazyk	Návod na použitie	PTCA dilatačný katéter	64
tr	Türkçe	Kullanma Talimatı	PTKA Dilatasyon Kateteri	66
vi	Tiếng Việt	Hướng dẫn sử dụng	Bóng nong PTCA	68

en Explanation of symbols on product labelling **de** Bedeutung der Symbole auf der Produktverpackung **fr** Explication des symboles sur l'étiquetage du produit **bg** Обяснение на символите върху етикетите на продуктите **zh** 产品标签上的符号说明 **cs** Význam symbolů na výrobku **da** Forklaring af symboler på produktmærkater **et** Toote siltidel esinevate sümbolite selgitused **es** Explicación de los símbolos en la ficha técnica del product **fi** Tuotteen pakkausmerkintöjen symbolien selitykset **el** Επεξήγηση των συμβόλων στην επισήμανση του προϊόντος **hr** Objašnjenje simbola na oznakama proizvoda **hu** A termék címkéjén megtalálható szimbólumok magyarázata **it** Legenda dei simboli sulle etichette dei prodotti **ko** 제품 라벨의 기호 설명 **lt** Gaminio etiketėje naudojamų simbolių paaiškinimas **lv** I zstrādājumu marķējumā iekļauto simbolu skaidrojums **mt** Spjegazzjoni tas-simboli fuq it-tikketta tal-prodott **nl** Verklaring van symbolen op productetiketten **no** Forklaring av symboler på produktmerkingen **pl** Objasnienie symboli na etykietcie produktu **pt** Explicação dos símbolos nos rótulos do produto **ro** Explicarea simbolurilor de pe eticheta produsului **ru** Объяснение символов на маркировке изделий **sr** Objašnjenje simbola na nalepnici proizvoda **sv** Förklaring av symboler på produktmärkningen **sl** Razlaga simbolov na nalepki izdelka **sk** Vysvetlenie symbolov na označení výrobku **tr** Ürün etiketlerindeki sembollerin açıklaması **vi** Giải thích các biểu tượng trên nhãn sản phẩm

 <p>Manufacturer / Hersteller / Fabricant / Производител / 制造商 / Výrobce / Producent / Tootja / Fabricante / Valmistaja / Κατασκευαστής / Proizvođač / Gyártó / Produttore / 제조사 / Gamintojas / Ražotājs / Manifattur / Fabrikant / Produzent / Producent / Fabricante / Fabricantul / Производител / Proizvođač / Tillverkare / Proizvajalec / Výrobca, Üretici / Nhà sản xuất</p>	<p>EC REP European authorized representative / Bevollmächtigter in der EU / Représentant agréé pour l'Europe / Европейски оторизиран представител / 欧洲授权代表 / Autorizovaný zástupce pro Evropu / Europæisk autoriseret repræsentant / Volitatud esindaja Euroopas / Representante europeo autorizado / Eurooppalainen valtuutettu edustaja / Εξουσιοδοτημένος αντιπρόσωπος για την Ευρώπη / Ovlašteni predstavnik za Evropu / Európai meghatalmazott képviselő / Rappresentante autorizzato per l'Europa / 유럽 공인 대리점 / Europos įgaliotasis atstovas / Pilnvarotais pārstāvis Eiropā / Rapprezentant Awtorizzat Ewropew / Gemachtigde vertegenwoordiger voor de Europese Unie / autoriseret representant i Europa / Upoważniony przedstawiciel w Unii Europejskiej / Representante autorizados na Europa / Reprezentantă autorizată în Europa / Уполномоченный представитель в Европе / Ovlašćeni zastupnik u Evropi / Europeisk auktoriserad representant / Evropski pooblašćeni zastopnik / Európsky autorizovaný zástupca / Avrupa Yetkili Temsilcisi / Đại diện được ủy quyền tại châu Âu</p>
 <p>EU Importer / EU-Importeur / Importateur pour l'UE / Вноситель от ЕС / 欧盟进口商 / Dovožce do EU / EU-importør / EL-i importtija / Importador de la UE / EU-tuoja / Εισαγωγέας για την ΕΕ / Uvoznik iz EU-a / EU-s behozatal / Importatore UE / EU 수입사 / ES importuotojas / ES importētājs / Importatur tal-UE / EU-Importeur / EU-importør / Importer do Unii Europejskiej / Importador na UE / Importator în UE / Импортёр в ЕС / Uvoznik za Evropsku uniju / EU-importör / Evropski uvoznik / Dovožca do EÚ / AB Ithalatçısı / Nhà nhập khẩu châu Âu</p>	<p>REF Catalogue number / Bestellnummer / Référence / Каталоген номер / 目录号 / Katalogové číslo / Katalognummer / Kataloogi number / Número de catálogo / Luettelonumero / Αριθμός καταλόγου / Kataloški broj / Katalógusszám / Codice a catalogo / 카탈로그 번호 / Katalogo numeris / Numurs katalogā / Numri tal-katalogu / Catalogusnummer / Katalognummer / Numer katalogowy / Referència / Număr de catalog / Номер по каталогу / Broj serije / Katalognummer / Kataloška številka / Katalogové číslo / Katalog numarası / Số danh mục</p>
<p>LOT Batch code / Chargenbezeichnung / N° de lot / Партиден код / 批号 / Kód šarže / Batchkode / Partii kood / Código de lote / Eränumero / Κωδικός παρτίδας / Šifra serije / Tételkód / Codice di lotto / 배치 코드 / Gaminio partijos kodas / Partijas kods / Kodici tal-lott / Partijnummer / Partikode / Kod partii / Código do lote / Codul seriei de fabricație / Код партии / Šifra partije / Batchkod / Koda serije / Číslo šarže / Parti kodu / Mã lô</p>	<p>UDI Unique Device Identifier / Eindeutige Produktidentifizierung / Identifiant unique du dispositif / Уникален идентификатор на изделие / 医疗器械唯一标识 / Jedinečný identifikátor zdravotníckého prostredku / Unik udstyrsidentifikation / Seadme kordumatu identifitseerimistunnus / Identificador único del dispositivo / Yksilöllinen laitetunniste / Μοναδικό αναγνωριστικό συσκευής / Jedinstvena identifikacija proizvoda / Egyedi eszközazonosító / Identificatore unico del dispositivo / 고유 장치 식별자 / Unikalus prietaiso identifikatorius / Unikālais ierīces identifikators / Identifikatur Uniku tat-Taġmir / Unieke identificatiecode van het hulpmiddel / Unik enhetsidentifikator / Niepowtarzalny identyfikator urządzenia / Identificador único do dispositivo / Identificator unic de dispozitiv / Уникальный идентификатор изделия / Jedinstveni identifikator uređaja / Unik produktidentifizierung / Eđinstveni identifikator pripomočka / Jedinečný identifikátor pomôcky / Benzersiz Cihaz Tanımlayıcı / Mã số nhận dạng thiết bị duy nhất</p>
<p>GTIN Global Trade Item Number / Global Trade Item Number / Numéro d'article commercial international / Глобален търговски номер на продукта / 全球交易品项识别代码 / Globální číslo obchodní položky / Globalt handelsvarenummer / Globaalne kaubaartikli number / Global Trade Item Number / Maailmanlaajuinen tuotenumero / Διεθνής κωδικός μονάδων εμπορίας / Globalni broj trgovačke jedinice / GTIN szám / Codice GTIN / 국제 거래 단위 번호 / Visuotiniai Prekių Numeriai / Globālais Tirdzniecības Vienības Numurs / Global Trade Item Number / Numru tal-Kummerċ Globali tal-Ogġett / Global Trade Item Number / Globalny numer jednostki handlowej / Número de referência de comércio global / Număr articol comercializat global / Глобальный номер предмета торговли / globalni broj trgovinske jedinice / GS1-artikelnummer Številka globalne trgovske enote / Globálne jedinečné číslo obchodnej položky / Global Ticari Ürün Numarası / Mã số phân định thương phẩm toàn cầu</p>	<p>Use-by date / Verwendbar bis / Date limite d'utilisation / Срок на годност / 使用截止日期 / Datum použitelnosti / Udløbsdato / Aegumiskuupeäev / Fecha de caducidad / Viimeinen käyttöpäivämäärä / Ημερομηνία λήξης / Rok upotrebe / Felhasználható / Data di scadenza / 유효일 / Data „Sunaudoti iki“ / Derīguma termiņš Data ta' skadenza / Uiterste gebruiksdatum / Siste forbruksdag / Data ważności / Prazo de validade / Data de expirare / Срок годности / Datum roka trajanja / Utgångsdatum / Rok uporabnosti / Použite do / Son kulanma tarihi / Hạn sử dụng</p>
 <p>Date of Manufacture / Herstellungsdatum / Date de fabrication / Дата на производство / 生产日期 / Datum výroby / Produktionsdato / Tootmiskuupeäev / Fecha de fabricación / Valmistuspäivämäärä / Ημερομηνία κατασκευής / Datum proizvodnje / A gyártás dátuma / Data di produzione / 제조일자 / Pagaminimo data / Izgatavošanas datums / Data ta' Produzzjoni / Fabricagedatum / Produksjonsdato / Data produkcji / Data de fabrico / Data fabricației / Дата изготовления / Datum proizvodnje /</p>	<p>Use-by date / Verwendbar bis / Date limite d'utilisation / Срок на годност / 使用截止日期 / Datum použitelnosti / Udløbsdato / Aegumiskuupeäev / Fecha de caducidad / Viimeinen käyttöpäivämäärä / Ημερομηνία λήξης / Rok upotrebe / Felhasználható / Data di scadenza / 유효일 / Data „Sunaudoti iki“ / Derīguma termiņš Data ta' skadenza / Uiterste gebruiksdatum / Siste forbruksdag / Data ważności / Prazo de validade / Data de expirare / Срок годности / Datum roka trajanja / Utgångsdatum / Rok uporabnosti / Použite do / Son kulanma tarihi / Hạn sử dụng</p>

Tillverkningsdatum / Datum izdelave / Dátum výroby /
Üretim Tarihi / Ngày sản xuất

STERILE EO

Sterilized using ethylene oxide / Sterilisiert mit Ethylenoxid / Stérilisé à l'oxyde d'éthylène / Стерилизован с използване на етиленов окис / 已使用环氧乙烷灭菌 / Sterilizováno etylénoxidem / Steriliseret med ethylenoxid / Steriliseeritud etüleenoksiidiga / Esterilizado con óxido de etileno / Steriloitu eteenioksidilla / Έχει αποστειρωθεί με χρήση οξειδίου του αιθυλενίου / Sterilizirano etilenoksidom / Etilén-oxidál sterilizálva / Sterilizzato a ossido di etilene / 산화에틸렌으로 멸균 처리됨 / Sterilizuota etileno oksidu / Sterilizēts ar etilēnoksidu / Sterilizat bl-užu ta' ethylene oxide / Gesteriliseerd met ethyleenoxide / Steriliseret ved hjelp av etylenoksid / Wysterylizowano tlenkiem etylenu / Esterilizado com óxido de etileno / Sterilizat cu etilenoxid / Стерилизовано окисью этилена / Sterilisanو pomoću etilen-okside / Steriliserad med etylenoxid / Sterilizirano z etilen oksidom / Sterilizované etylénoxidom / Etilen oksitle sterilize edilmıştır / Tiệt trùng bằng etilen oxit

Do not use if package is damaged / Bei beschädigter Verpackung nicht verwenden / Ne pas utiliser si l'emballage est endommagé / Да не се използва, ако опаковката е повредена / 如包装损坏请勿使用 / Je-li obal poškozený, výrobek nepoužívejte / Må ikke anvendes, hvis pakningen er beskadiget / Mitte kasutada, kui pakend on kahjustatud / No utilizar si el envase está dañado / Ei saa kasutada, jos pakkaus on vahingoittunut / Μη χρησιμοποιείτε σε περίπτωση που η συσκευασία έχει υποστεί ζημιά / Ne upotrebljavati ako je pakovanje oštećeno / Ne használja fel, ha a csomagolás sérült / Non utilizzare se la confezione è danneggiata / 포장이 손상된 경우 사용 금지 / Nenaudoti, jei pakuotė pažeista / Nelietot, ja iepakojums ir bojāts / Tužax jekk il-pakkett ikun fih il-ħsara / Niet gebruiken als de verpakking beschadigd is / Skal ikke brukes hvis emballasjen er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Não utilizar se a embalagem estiver danificada / Nu se utiliza dacă ambalajul este deteriorat / Не използуйте в случае повреждения упаковки / Не користити ако је амбалажа оштећена / Använd inte om förpackningen är skadad / Ne uporabljajte, če je embalaža poškodovana / Nepoužívajte, ak je obal poškodený / Ambalaj hasarlıysa kullanmayın / Không sử dụng nếu bao bì bị hỏng

Single sterile barrier system with protective packaging outside / Einfaches Sterilbarriersystem mit äußerer Schutzverpackung / Système de barrière stérile simple avec emballage protectif extérieur / Единична стерилна бариерна система със защитна опаковка отвън / 带保护性外包装的单一无菌屏障系统 / Systém jedné sterilní bariéry s vnějším ochranným balením / System med enkelt sterilbarriere med beskyttende udvendig emballage / Ühekordne steriilne barjäärisüsteem kaitsva välispakendiga / Sistema de barrera estéril único con embalaje protector exterior / Steriilin esteen järjestelmä ja suojaava ulkorakkaus / Σύστημα μονού στέρπου φραγμού με προστατευτική συσκευασία στο εξωτερικό / Sustav jednostruke sterilne barijere sa zaštitnim vanjskim pakiranjem / Egyszeres sterilgátrendszer külső védőcsomagolással / Sistema di barriera sterile singola con imballaggio protettivo esterno / 외부에 보호 포장에 있는 단일 멸균 장벽 시스템 / Vieno sterilaus barjero sistema su apsaugine išorine pakuote / Vienas sterila barjeras sistēma ar aizsargiekpakojumu ārpusē / Sistema sterili ta' barriera unika b'imballaġġ protettiv fuq barra / Systeem met enkele steriele barrière met beschermende verpakking aan de buitenkant / Enkelt sterilt barriersystem med beskyttende innpakning / System pojedynczej bariery sterylnej z osłoną ochronną na zewnątrz / Sistema de barreira estéril simples com embalagem exterior de proteção / Sistem de barieră sterilă unică cu ambalaj de protecție în exterior / Одинарная стерильная барьерная система со внешней защитной упаковкой / Jednostruki sterilni sistem barijere sa zaštitnim spoljnim pakovanjem / Enkelt sterilbarriärsystem med skyddsförpackning utanför / Sistem enojne sterilne pregrade z zaščitno zunanjo embalažo / Jednoduchý systém sterilnej bariéry s vonkajším ochranným obalom / Diš koruyucu ambalajlı tek steril bariyer sistemi / Hệ thống bảo vệ vô khuẩn một lớp có màng bảo vệ bên ngoài

One device in a box / Ein Produkt im Karton / Un dispositif dans une boîte / Одно изделие в кутия / 每盒一个设备 / Jeden kus v krabici / En enhed i en æske / Karbis on üks seade / Un dispositivo en una caja / Yksi laite pakkauksessa / Μία συσκευή σε ένα κουτί / Jedan proizvod u kutiji / Egy eszköz van a dobozban / Un dispositivo per confezione / 상자 하나에 기구 한 개 포함 / Déžutéje yra vienas prietaisas / Viena ierīce kārbā / aghmir wiehed f'kaxxa / Eén hulpmiddel per doos / Én enhet i en eske / Jedno urządzenie w opakowaniu / Um dispositivo por caixa / Un dispozitiv într-o cutie / В упаковке одно изделие / Jedan uređaj u kutiji / En enhet i en ask / En pripomoček v škatli / Jedno zariadenie v krabici / Kutuda bir Cihaz / Một thiết bị trong một hộp

Consult instructions for use / Gebrauchsanweisung beachten / Consulter le mode d'emploi / Консультируйтесь с инструкцией за употреба / 请参阅使用说明 / Přečtěte si návod k použití / Se brugsanvisningen / Lugege kasutusjuhendit / Consulte las instrucciones de uso / Lue käyttöohjeet / Συμβουλευτείτε τις οδηγίες χρήσης / Pročitajte upute za upotrebu / Olvassa el a Használati útmutatót / Consultare le istruzioni per l'uso / 사용 설명서 참조 / Susipažinkite su naudojimo instrukcijomis / Skatīt lietošanas pamācību / Ikkonsultala l-istruzzjonijiet għall-użu / Raadpleeg de gebruiksaanwijzing / Se bruksanvisningen / Zapoznać się z instrukcją obsługi / Consultar as instruções de utilização / Consultați instrucțiunile de utilizare / Обратитесь к инструкции по применению / Pogledati uputstvo za upotrebu / Se bruksanvisningen / Glejte navodila za uporabo / Prečítajte si návod na použitie / Kullanna talimatna başvurun / Tham khảo hướng dẫn sử dụng

Caution / Vorsicht / Avertissement / Внимание / 警告 / Pozor / Forsigtig / Ettevaatust / Precaución / Huomio / Προσοχή / Oprez / Figyelem! / Attenzione / 주의 / Atsargiai / Uzmanību! / Attenzjoni / Voorzichtig / Forsiktig / Przestroga / Cuidado / Atenție / Внимание! / Oprez / Var försiktig / Pozor / Upozornenie / Dikkat / Thận trọng

Do not re-use / Nicht zur Wiederverwendung / Ne pas reutiliser / Да не се използва повторно / 请勿重复使用 / Nepoužívat opakovaně / Må ikke genbruges / Mitte kasutada korduvalt / No reutilizar / Ei saa kasutada uudelleen / Μη επαναχρησιμοποιείτε / Ne upotrebljavati ponovno / Ne használja fel újra /


Do not resterilize / Nicht erneut sterilisieren / Ne pas restériliser / Да не се стерилизира повторно / 请勿重新灭菌 / Nesterilizovat opakovaně / Må ikke resteriliseres / Mitte steriliseerida uuesti / No reesterilizar / Ei saa steriloida uudelleen / Μη επαναποστειρώνετε / Ne sterilizirati ponovno / Ne

Non riutilizzare / 재사용 금지 / Nenaudoti pakartotinai / Nelietot atkārtoti / Tužax mill-ġdid / Niet opnieuw gebruiken / Skal ikke gjenbrukes / Nie używać ponownie / Não reutilizar / A nu se reutiliza / Только для одноразового применения / Ne koristiti ponovo / Återanvänd inte / Ne uporabljajte ponovno / Nepoužívejte opakovane / Tekrar kullanmayın / Không tái sử dụng

sterilizálja újra! / Non risterilizzare / 재멸균 금지 / Nesterilizuoti pakartotinai / Nesterilizēt / Tisterilizax mill-ġdid / Niet opnieuw steriliseren / Skal ikke steriliseres på nytt / Nie sterylizować ponownie / Não voltar a esterilizar / A nu se reesteriliza / Повторная стерилизация запрещена / Ne sterilisasi ponovo / Får inte omsteriliseras / Ne sterilizirajte ponovno / Nesterilizujte opakovane / Tekrar sterilize etmeyin / Không tiệt trùng lại

Keep away from sunlight / Vor Sonnenlicht schützen / Protéger du rayonnement solaire / Да се пази от слънчева светлина / 避免日光照射 / Uchovávejte mimo dosah slunečního záření / Må ikke udsættes for sollys / Hoida otsese päikesevalguse eest / Evitar la exposición a la luz del sol / Suojattava auringonvaloita / Φυλάσσετε μακριά από το ηλιακό φως / Držati podalje od sunčeve svjetlosti / Napfénytől védve tárolandó / Tenere al riparo dalla luce solare / 직사광선을 피할 것 / Saugoti nuo saulės spindulių / Neturēt saules gaismā / Zomm 'il bogħod mid-dawl tax-xemx / Beschermen tegen zonlicht / Må holdes unna sollys / Chronić przed światłem słonecznym / Manter afastado da luz solar / A se feri de lumina solară / Беречь от попадания солнечных лучей / Držati dalje od izvora svetlosti / Förvaras skyddat från solljus / Hranite stran od sončne svetlobe / Chránite pred slnečným žiarením / Güneş ışığından uzak tutun / Tránh xa ánh sáng mặt trời

Store between 15 and 30 °C / Bei 15 bis 30 °C lagern / Conserver à une température comprise entre 15 et 30 °C / Да се съхранява при температура между 15 и 30 °C / 存放温度为 15 - 30°C / Uchovávat mezi 15 a 30 °C / Opbevares mellem 15 og 30 °C / Hoiustada temperatuuril vahemikus 15 kuni 30 °C / Conserver entre 15 y 30 °C / Säilytetään 15–30 °C:ssa / Φυλάσσετε μεταξύ 15 και 30 °C / Čuvati na temperaturi od 15 do 30 °C / 15 és 30 °C között tárolandó / Conservare tra 15 e 30 °C / 15 ~ 30°C, 에서 보관 / Laikyti esant temperatūrai nuo 15 iki 30 °C / Uzglabāt temperatūrā no 15 līdz 30 °C / Añzen f'temperatura bejn 15 u 30°C / Bewaren tussen 15 en 30 °C / Oppbevares ved mellom 15 og 30 °C / Przechowywać w temperaturze 15–30°C / Conserver entre 15 e 30 °C / A se păstra la temperaturi între 15 și 30°C / Хранить при температуре 15–30 °C / Čuvati na temperaturi između 15 i 30 °C / Förvaras i temperaturer mellan 15 och 30 °C / Hranite pri temperaturi med 15 in 30 °C / Skladujte pri teplote medzi 15 a 30 °C / 15 ile 30°C arasında saklayın / Bảo quản ở nhiệt độ từ 15 đến 30°C

Keep dry / Trocken aufbewahren / Conserver à l'abri de l'humidité / Да се пази в сухо състояние / 保持干燥 / Uchovávejte v suchu / Opbevares tørt / Hoida kuivana / Conservar en un lugar seco / Säilytetävä kuivassa / Διατηρείτε στεγνό / Čuvati na suhom / Szárazon tartandó / Conservare in luogo asciutto / 건조한 곳에 보관 / Laikyti sausai / Sargāt no mitruma / Zomm il-prodott xott / Droog bewaren / Må holdes tørr / Przechowywać w suchym miejscu / Manter seco / A se păstra uscat / Хранить в сухом месте / Držati na suvom mestu / Förvaras tort / Hranite na suhem / Udržujte v suchu / Kuru tutun / Giữ khô ráo

Medical device / Medizinprodukt / Dispositif médical / Медицинско изделие / 医疗器械 / Zdravotnický prostředek / Medicinsk udstyr / Meditsiiniseade / Dispositivo medico / Lääkinnällinen laite / Ιατρο-τεχνολογικό προϊόν / Medicinski proizvod / Orvos-technikai eszköz / Dispositivo medico / 의료 장치 / Medicinos prietaisai / Medicīniska ierīce / Tagħmir mediku / Medisch hulpmiddel / Medisinsk utstyr / Urządzenie medyczne / Dispositivo medico / Dispozitiv medical / Медицинское изделие / Medicinski uređaj / Unik produktidentifizierung / Medicinski pripomoček / Zdravotnícka pomôcka / Tibbi cihaz / Thiết bị y tế

Balloon Length / Ballonlänge / Longueur du ballonnet / Дължина на балона / 球囊长度 / Délka balónku / Ballonlængde / Ballooni pikkus / Longitud del balón / Pallon pituus / Μήκος μπαλονιού / Duljina balona / Ballon hosszúsága / Lunghezza del palloncino / 풍선 길이 / Balionélio ilgis / Balona garums / Tul tal-Buzzieqa / Ballonlengte / Lengde på ballong / Długość balonu / Comprimento do balão / Lungimea balonului / Длина баллона / Dužina balona / Ballonglängd / Dolžina balona / Dĺžka balónika / Balon Uzunluğu / Chiều dài bóng

Balloon Diameter / Ballondurchmesser / Diamètre du ballonnet / Диаметър на балона / 球囊直径 / Průměr balónku / Ballondiameter / Ballooni läbimõõt / Diámetro del balón / Pallon halkaisija / Διαμέτρος μπαλονιού / Promjer balona / Ballon átmérője / Diametro del palloncino / 풍선 직경 / Balionélio skersmuo / Balona diametrs / Dijametru tal-Buzzieqa / Ballondiameter / Diameter på ballong / Šrednica balonu / Diâmetro do balão / Diametru balonului / Диаметр баллона / Prečnik balona / Ballongdiameter / Premer balona / Priemer balónika / Balon Çapı / Đường kính bóng

Recommended Guide Wire / Empfohlener Führungsdraht / Fil de guidage recommandé / Препоръчан водач / 推荐导丝 / Doručovaný vodící drát / Anbefalet ledetråd / Soovitav juhtraat / Alambre guía recomendado / Suositeltu ohjausvaijeri / Συνιστώμενο οδηγό σύρμα / Preporučena žica vodilica / Ajánlott vezetődó / Filo guida consigliato / 권장 가이드 와이어 / Rekomenduojama kreipiamoji viela / Ieteicamā vadītājstīga / Wajer Gwida Rakkomandat / Aanbevolen voerdraad / Anbefalt guidewire / Zalecana średnica prowadnika / Fio guida recomendado / Fir de ghidaj recomandat / Рекомендуемый проволочный проводник катетера / Preporučena žica-vodica / Rekomenderad ledare / Priporučena vodilna žica / Odporučaný vodiací drôt / Önerilen Kılavuz Tel / Dây dẫn được khuyến nghị

Minimum Guiding Catheter / Mindestgröße Führungskatheter / Cathéter de guidage minimal / Минимален водещ катетър / 最小导引导管 / Minimální zaváděcí katétr / Minimumsledekateter / Minimaalne juhtkateeter / Catéter guía mínimo / Ohjauskatetrin vähimmäishalkaisija / Ελάχιστο μέγεθος οδηγού καθετήρα / Minimalni vodeći kateter / Minimalis vezetőkátéter / Catetere guida minimo / 최소 가이드팅 카테터 / Mažiausias kreipiamasis kateteris / Minimālais vadītājkatētrs / Kateter Gwida Minim / Minimale geleidekatheter / Minste ledekateter / Minimalna średnica cewnika prowadzącego / Cateter guida mínimo / Cateter de ghidaj minim / Минимальный направляющий катетер / Minimalni uvodni kateter / Minsta styrkateter / Najmanjši vodilni kateter / Minimálny vodiaci katéter / Minimum Kılavuz Kateter / Ống thông dẫn hướng có kích thước tối thiểu

Inflation Pressure / Inflationsdruck / Pression d'insufflation / Налягане на раздуване / 扩张压力 / Plnicí tlak / Inflationstryk / Taitmisröhk / Presión de

Nominal Pressure / Nenndruck / Pression nominale / Номинално налягане / 标称压力 / Jmenovitý tlak / Nominel tryk / Nimiröhk / Presión nominal /

<p>inflado / Τάυττόρρarine / Πίεση πλήρωσης / Tlak naruhiwanja / Dilatációs nyomás / Pressione di gonfiaggio / 팽창압력 / Išpūtimo slėgis / Piepildīšanas spiediēns / Pressjoni tal-Inflazzjoni / Vuldruk / Oppblåsningstrykk / Ciśnienie napełniania / Pressão de enchimento / Presiune de umflare / Давление расширения / Pritisak naduvavanja / Uppblåsningstryck / Tlak narihovanja / Tlak nafúknútia / Şişirme Basıncı / Áp suất bơm phồng</p>	<p>Nimellispaine / Ονομαστική πίεση / Nazivni tlak / Névleges nyomás / Pressione nominale / 공칭압력 / Nominalusis slėgis / Nominālais spiediens / Pressjoni Nominali / Nominale druk / Nominelt trykk / Ciśnienie nominalne / Pressão nominal / Presiune nominală / Номинальное давление / Nominalni pritisak / Nominellt tryck / Nominálny tlak / Nazivni tlak / Nominal Basınç / Áp suất danh định</p>
<p>RBP Rated Burst Pressure / Nennberstdruck / Pression d'éclatement nominale / Номинално налягане на спукване / 额定爆破压力 / Jmenovitý destrukční tlak / Nominelt sprængningstryk / Hinnanguline lõhkemisrõhk / Presión nominal de estallido / Nimellinen murtumisaine / Ονομαστική πίεση ρήξης / Nazivni tlak rasprskavanja / Névleges szétrepedési nyomás / Pressione nominale di scoppio / 정격파열압력 / Apskaičiuotasis išpūtimo slėgis / Nominālais pārraušanas spiediens / Pressjoni Nominali tat-Tifqigh / Nominale barstdruk / Nominelt bursttrykk / Nominalne ciśnienie rozerwania / Pressão de rebentamento nominal / Presiune nominală de spargere / Расчетное давление разрыва / Označen pritisak pucanja balona / Nominellt bristningstryck / Nazivni porušitveni tlak / Menovitý tlak pretrhnútia / Anma Patlama Basıncı / Áp suất vỡ bóng định mức</p>	<p>DEFLATION Deflation Time / Deflationsdauer / Temps de dégonflement / Време за отдуване / 收缩时间 / Doba vypouštění / Deflationstid / Tühjendamisaeg / Tiempo de desinflado / Τυήjenemisaika / Χρόνος σμύπτωσης / Vrijeme ispuhivanja / Leengedési idő / Tempo di sgonfiaggio / 수축 시간 / Išleidimo laikas / Gaisa izlaišanas ilgums / Īin ta' deflazzjoni / Ledigingstijd / Tømmetid / Czas utraty powietrza / Tempo de desinsuflação / Timp de deflație / Время сдувания / Vreme izduvananja balona / Tömningstid / Čas deflacije / Čas vyfúknútia / Īnme Süresi / Thời gian xả xẹp</p>
<p> hydrophilic Hydrophilic Coating / Hydrophile Beschichtung / Revêtement hydrophile / Хидрофилно покритие / 亲水涂层 / Hydrofilní povlak / Hydrofil belægning / Hüdروفилне кате / Revestimiento hidrófilo / Hydrofilinen pinnoite / Υδρόφιλη επικάλυψη / Hidrofilni premaz / Hidrofil bevonat / Rivestimento idrofilico / 친수 코팅 / Hidrofiliné danga / Hidrofilis pārklājums / Kisi Idrofiliku / Hydrofiele coating / Hydrofilt belegg / Powłoka hydrofilowa / Revestimento de hidrocoloide / Īnveliř hidrofil / Гидрофильное покрытие / Hidrofilni premaz / Hydrofil belægning / Hidrofilna prevleka / Hydrofilný povlak / Hidrofilik Kaplama / Lóp phủ hút nước</p>	<p> hydrophobic Hydrophobic Coating / Hydrophobe Beschichtung / Revêtement hydrophobe / Хидрофобно покритие / 疏水涂层 / Hydrofobní povrchová úprava / Vandskyende belægning / Hüdروفобне каттекиht / Revestimiento hidrofóbico / Vettähylkivä pinnoite / Υδρόφοβη επικάλυψη / Hidrofobni premaz / Hidrofób bevonat / Rivestimento idrofobico / 친수 코팅 / Hidrofobinē danga / Hidrofobs pārklājums / Kisi Idrofobiku / Hydrofobe coating / Vannavvisende belegg / Powłoka hydrofobowa / Revestimento hidrófugo / Īnveliř hidrofob / Гидрофобное покрытие / Hidrofobni premaz / Hydrofob belægning / Hidrofobna prevleka / Hydrofóbny povlak / Hidrofobik Kaplama / Lóp phủ không thấm nước</p>
<p> 1 One proximal marker / ein proximal Marker / un gerèpe proximal / един проксимален маркер / 一个近端标记 / jeden proximální marker / en proksimal markør / üks proksimaalne tähis / un marcador proximal / yksi proksimaalinen merkki / ένας εγγύς δείκτης / jedan proksimalni marker / egy proximális marker / un marker prossimale / 근위 마커 한 개 / viena žyma proksimaliojoje dalyje / Viens proksimālais marķieris / markatur prossimali wiehed / één proximale markering / én proksimal markør / jeden znacznik proksymalny / um marcador proximal / un marker proximal / один проксимальный маркер / jedan proksimalan marker / en proximal markør / ena proksimalna oznaka / jedna proximálna značka / bir proksimal işaret / đầu một ở đầu gần</p>	<p> 1 One middle marker / ein zentraler Marker / un repère central / един среден маркер / 一个中间标记 / jeden středový marker / en mellemmarkør / üks keskmine tähis / un marcador intermedio / yksi keskikohdan merkki / ένας μεσσίος δείκτης / jedan srednji marker / Egy középső marker / un marker intermedio / 중간 마커 한 개 / viena žyma vidurinēje dalyje / viens vidējais marķieris / markatur tan-nofs wiehed / één middenmarkering / én midtre markør / jeden znacznik środkowy / um marcador medial / un marker median / один средний маркер / jedan midalni marker / en mittmarkør / ena sredinska oznaka / jedna stredová značka / bir orta işaret / đầu một ở giữa</p>
<p>Rx only USA only: Caution: Federal law restricts the device for sale by or on the order of a physician</p>	
<p>NOTE / HINWEIS / REMARQUE / ЗАБЕЛЕЖКА / 注意 / POZN / BEMÆRK / MÄRKUS / NOTA / HUOMAA / ΣΗΜΕΙΩΣΗ / NAPOMENA / MEGJEGYZÉS / NOTA / 참고 / PASTAVA / PIEZĪME / NOTA / NB / MERK / UWAGA / NOTA / NOTĀ / ПРИМЕЧАНИЕ / NAPOMENA / OBS! / ΟΡΟΜΒΑ / POZNÁMKA / NOT / LU'U Ý</p>	<p>1 atm ≅ 100 kPa = 1 Bar</p>

Compliance Chart / Compliance-Tabelle / Tableau de compliance / Таблица за съответствие / 顺应性表 / Tabulka compatibility / Overensstemmingschema / Ühilduvuse tabel / Tabla de compliance / Vaatimustenmukaisuus-listaus / Πίνακας συμμόρφωσης / Tablica skladnosti / Megfelelési diagram / Tabella di conformità / 규정 차트 / Atitikčių lentelė / Atbilstības tabula / Čart tal-Konformità / Compliantietabel / Samsvarsdiagram / Karta podatności / Quadro de distensibilidade / Grafic de conformitate / Таблица соответствия / Tablica usklađenosti / Tabell över överensstämmelse / Tabela skladnosti / Tabuľka poddajnosti / Esneklik Tablosu / Bảng độ giãn nở

	(IP) _[atm]	Ø Balloon [mm]							(IP) _[atm]
		1.5	2.0	2.5	3.0	3.5	4.0	4.5	
OPN NC	10 NOM	1.50	2.00	2.50	3.00	3.50	4.00	4.50	10 NOM
	20	1.55	2.10	2.60	3.14	3.67	4.19	4.66	20
	30	1.60	2.18	2.70	3.29	3.85	4.37	4.82	30
	35 RBP	1.68	2.22	2.77	3.36	3.91	4.44	4.90	35 RBP
BEO NC	10 NOM	1.50	2.00	2.50	3.00	3.50	4.00	4.50	10 NOM
	12	1.56	2.06	2.54	3.04	3.57	4.10	4.60	12
	14	1.58	2.10	2.57	3.07	3.61	4.17	4.67	14
	16	1.59	2.11	2.60	3.11	3.65	4.22	4.71	16
	18	1.62	2.13	2.63	3.17	3.69	4.27	4.77	18
	20	1.64	2.16	2.67	3.22	3.73	4.31	4.83	20
	21 RBP	1.65	2.18	2.69	3.24	3.75	4.34	4.85	21 RBP
EasyT	10 NOM	1.55	2.02	2.45					
	12	1.59	2.11	2.54	3.01	3.43	3.95	4.45	12 NOM
	14	1.63	2.17	2.61	3.12	3.55	4.11	4.61	14
	16	1.66	2.21	2.66	3.20	3.65	4.23	4.73	16
	18	1.68	2.24	2.71	3.26	3.72	4.32	4.82	18
	20	1.70	2.27	2.76	3.31	3.79	4.39	4.89	20
	21 RBP	1.73	2.29	2.78	3.35	3.83	4.44	4.93	21 RBP
NIC Nano hydro		0.85							
	10 NOM	0.85							
	12	0.88							
	14	0.90							
	16	0.92							
	18	0.94							
	20	0.96							
	21 RBP	0.97							
NIC 1.1 hydro		1.1							
	10 NOM	1.10							
	12	1.14							
	14	1.18							
	16	1.22							
	18	1.26							
	20	1.30							
21 RBP	1.34								

Do not exceed Rated Burst Pressure / Nennberstdruck nicht überschreiten / Ne pas dépasser la pression d'éclatement nominale / Не превишавайте номиналното налягане на спукване / 请勿超过额定爆破压力 / Nepřekračujte nominální destrukční tlak / Den nominelle øvre trykgrænse må ikke overstiges / Mitte ületada hinnangulist lõhkemisrõhku / No sobrepasar la presión nominal de estallido / Άλἄ ὑλιτἄ nimellistἄ murtumispainetta / Μην υπερβαίνετε την ονομαστική πίεση ρήξης / Ne premašivati nazivni tlak rasprskavanja / Ne lépje túl a névleges szétrepedési nyomást / Non superare la pressione nominale di scoppio / 정격파열압력을 초과하지 마십시오 / Neviršykite nominalaus plyšimo slėgio / Nepārsniedziet nominālo pārraušanas spiedienu / M'għandekx taqbeż il-Pressjoni ta' Fqigh Gradat / Nominale barstdruk niet overschrijden / Nominelt bursttrykk må ikke overskrides / Nie przekraczać nominalnego ciśnienia rozerwania / Não ultrapassar a pressão nominal de rebentamento / Nu depășiți presiunea nominală de spargere / Не превышайте расчетное давление разрыва / Nemojte da prekoračite označen pritisak pucanja balona / Överskrid inte det nominella bristningstrycket / Ne presežite nazivnega porušitvenega tlaka / Neprekrāčujtie menovitý tlak pretrhnutia / Anma Patlama Basıncını geçmeyin / Không vượt quá Áp suất vỡ bóng định mức

DESCRIPTION

The percutaneous transluminal coronary angioplasty (PTCA) Dilatation Catheter operates on the principle of hydraulic pressurization applied through an inflatable balloon attached to the distal end. Distally it has a double lumen shaft housing the guide wire lumen and the inflation media lumen. The balloon is designed to uniformly expand to the nominal diameter and length at nominal pressure as labeled. The device may be expanded to pressures above the nominal pressure, but not above the rated burst pressure. The device features radiopaque markers (one or two depending on the size) to enable fluoroscopic visualization of the balloon during use. Furthermore, a brachial and a femoral marker are located on the shaft. The balloon catheter is to be used with guide wires with maximum 0.014" diameter in rapid exchange configuration. A hydrophobic coating is applied to the distal section of the catheter to increase lubricity and improve usability.

Caution: Federal Law (USA) restricts this device for sale by or on the order of a physician.

CONTENTS

One PTCA Dilatation Catheter, one instructions for use.

INDICATIONS FOR USE

The PTCA Dilatation Catheter is indicated for:

- Dilatation of coronary artery stenosis
- Dilatation of bypass graft stenosis
- Post dilatation of balloon expandable coronary stents

INTENDED USE

The PTCA Catheter is intended for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. The balloon dilatation catheter is also indicated for post deployment expansion of balloon expandable coronary stents.

TARGET USER GROUP

PTCA Dilatation Catheters should only be used by physicians who are experienced in coronary angiography and the PTCA technique.

TARGET PATIENT POPULATION

The product shall only be used in adult persons above 18 years of age.

CONTRAINDICATIONS

The PTCA Dilatation Catheter is contraindicated for:

- non-candidates for coronary bypass graft surgery
- lesions that cannot be reached or treated with the system
- unprotected left main coronary artery
- coronary artery spasm in the absence of significant stenosis

WARNINGS

- **STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.
- For single patient, single procedure use only. **Do not resterilize or reuse the device.** Resterilization or reuse may compromise device performance and increase the risk of inappropriate resterilization and cross contamination.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- For patients who are not acceptable candidates for coronary artery bypass graft surgery careful consideration is required, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis. Consider both balloon diameters if kissing balloon technique is performed.
- Do not use the PTCA Dilatation Catheter if its package has been opened or damaged.
- Do not use or straighten a PTCA Dilatation Catheter if the shaft is kinked as this may lead to breakage of the shaft.
- Do not expose the PTCA Dilatation Catheter to any organic solvents, e.g. alcohol as they may have an impact on the properties of the PTCA Dilatation Catheter.
- Use only appropriate balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.
- Do not advance or retract the PTCA Dilatation Catheter if unusual resistance is felt. When excessive force needs to be applied to move the PTCA Dilatation Catheter, stop the procedure before proceeding. Determine the root cause of the resistance before moving the catheter. Advancement by force may result in damage to the vessel and/or the PTCA Dilatation Catheter. This may necessitate recovery of fragments.
- Do not advance or retract the PTCA Dilatation Catheter unless the balloon is fully deflated under vacuum.
- Balloon pressure should not exceed the rated burst pressure (RBP) indicated on the package label for each balloon. The RBP is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over-pressurization. If inflated above the rated burst pressure, there is a risk of catheter burst.
- Do not torque the PTCA Dilatation Catheter more than one (1) full turn.
- Bench testing was conducted with a standard balloon expandable stent, e.g. Biotronik Prokinetic. Consideration should be taken when the device is used with different manufacturer's stents due to differences related to the stent design, inner surface and other characteristics. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.
- The balloon may be inflated several times, but it should not be inflated more than three times to the specified rated burst pressure. If inflated more than three times to the rated burst pressure, there is a risk of catheter burst.
- When the PTCA Dilatation Catheters is in the body, it should only be manipulated while under fluoroscopic observation with radiographic equipment that provides high quality images.

PRECAUTIONS

- PTCA Dilatation Catheters should only be used by physicians who are experienced in coronary angiography and the PTCA technique.
- During the procedure, appropriate anticoagulation and antiplatelet administration must be provided to the patient. Coronary vasodilator therapy is usually provided to the patient. Antiplatelet aggregation treatment and vasodilators are usually continued in all cases for a period of time determined by the physician.
- Prior to angioplasty, the catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.
- If resistance during insertion of the PTCA Dilatation Catheter is encountered do not use force to advance the catheter as this may result in damage to the catheter.
- If resistance during removal of the PTCA Dilatation Catheter is encountered do not use force to retract the catheter as this may result in damage to the catheter. Simultaneous removal of the entire system (e.g. guide wire and guiding catheter) might be advised.
- The safety and effectiveness of this PTCA Dilatation Catheter for the treatment of in-stent restenosis have not been established.
- Use the PTCA Dilatation Catheter prior to the "Use-by Date" specified on the package.

ADVERSE EFFECTS AND RISKS OF USE

These can include and are not be limited to:

- Acute myocardial infarction,
- Acute vessel closure,
- Arteriovenous fistula,
- Arrhythmias including ventricular fibrillation,
- Coronary artery spasm,
- Coronary vessel dissection, perforation, rupture or other injury,
- Death,
- Drug reaction or allergic reaction to contrast medium, antiplatelets, anticoagulants,
- Embolism,
- Hemorrhage or hematoma,
- Hypo- / hypertension,
- Infection,
- Restenosis of the dilated artery,
- Stroke, air embolism and embolism due to fragmentation of thrombotic or atherosclerotic material,
- Total occlusion of the coronary artery or bypass graft,
- Unstable angina.

DIRECTIONS FOR USE

PREPARATION OF THE PTCA DILATATION CATHETER

To properly prepare the catheter for use, the following steps should be completed:

1. Remove PTCA Dilatation Catheter gently out of packaging (sterile packaging and dispenser).
2. Fill a plastic syringe with 6 to 8 ml of inflation medium (e.g., 50:50 sterile mixture of a contrast medium and saline). Attach the syringe to the luer lock connector located on the proximal end of the PTCA Dilatation Catheter.
3. Remove the PROTECTIVE SHEATH from the balloon. Assure that PACKAGING WIRE is also removed.
NOTE: Do not use excessive force to remove the protective sheath or packaging wire.
4. Pull a vacuum on the syringe to remove any air in the catheter lumen or in the balloon. Gentle tapping of the PTCA Dilatation Catheter may help to remove air from the balloon and the PTCA Dilatation Catheter into the syringe.
5. Flush the guide wire lumen of the PTCA Dilatation Catheter with saline solution.

INSERTION TECHNIQUE

6. Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.
7. Advance a guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the PTCA Dilatation Catheter over the proximal end of the guide wire. Ensure that the guide wire exits the PTCA Dilatation Catheter through the guide wire exit port.
8. The hemostasis valve of the guiding catheter should be gradually tightened to control back flow.
NOTE: Excessive valve tightening may affect balloon inflation / deflation time as well as movement of the guide wire.
9. Track the PTCA Dilatation Catheter over the guide wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion under fluoroscopy.
PRECAUTION: If resistance during insertion is encountered do not use force to advance the PTCA Dilatation Catheter as this may result in damage to the PTCA Dilatation Catheter. Stop the procedure and determine the cause.
NOTE: Do not torque the PTCA Dilatation Catheter more than one (1) full turn.

BALLOON INFLATION

10. Inflate the balloon first to a low pressure to confirm location. Once proper balloon position is set under fluoroscopy, high pressure inflation may be initiated.
NOTE: Do not exceed the rated burst pressure indicated on the label.
11. After the first inflation and subsequent deflation, assess distal coronary blood flow by fluoroscopy (contrast medium injection through the guiding catheter while the balloon is deflated).
NOTE: A significant stenosis may require successive inflations. If necessary, the PTCA Dilatation Catheter may be exchanged without removing the guide wire and without an exchange wire.

CATHETER REMOVAL

12. After dilation, fully deflate the balloon and confirm that the balloon is fully deflated. Apply negative pressure until the balloon is fully deflated as observed under fluoroscopy.
13. Loosen the hemostatic valve on the guiding catheter. Hold the guide wire position and withdraw the PTCA Dilatation Catheter into the guiding catheter. Remove the PTCA Dilatation Catheter from the system.
PRECAUTION: If resistance is encountered during removal do not use force to retract the device as this may result in damage to the PTCA Dilatation Catheter. Simultaneous removal of the entire system (e.g. PTCA Dilatation Catheter, guide wire and guiding catheter) might be advised.
14. Perform a post angioplasty angiogram by contrast medium injection through the guiding catheter.
NOTE: Inspect the balloon and the PTCA Dilatation Catheter for any potential damage after wiping clean with gauze soaked with sterile normal saline.
15. Dispose of PTCA Dilatation Catheter as per internal hospital or national requirements and regulations. In case of any damage to the PTCA Dilatation Catheter of unknown cause, return the device to the manufacturer.

REFERENCES

The physician should consult recent literature on current medical practice on PTCA, such as published by the American College of Cardiology, the American Heart Association or European Society of Cardiology.

DISCLAIMER

ALTHOUGH THE PTCA DILATATION CATHETER HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, SIS MEDICAL AG AND ITS AFFILIATES HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. SIS MEDICAL AG AND ITS AFFILIATES, THEREFORE, DISCLAIM ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PTCA DILATATION CATHETER, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. SIS MEDICAL AG AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PTCA DILATATION CATHETER, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND SIS MEDICAL AG AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PTCA DILATATION CATHETER.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Any serious incident involving this product shall be reported to the manufacturer and to the competent authority of the country where the product was used.

DESCRIPTION

The percutaneous transluminal coronary angioplasty (PTCA) Dilatation Catheter operates on the principle of hydraulic pressurization applied through an inflatable balloon attached to the distal end. Distally it has a double lumen shaft housing the guide wire lumen and the inflation media lumen. The balloon is designed to uniformly expand to the nominal diameter and length at nominal pressure as labeled. The device may be expanded to pressures above the nominal pressure, but not above the rated burst pressure. The device features radiopaque markers (one or two depending on the size) to enable fluoroscopic visualization of the balloon during use. Furthermore, a brachial and a femoral marker are located on the shaft. The balloon catheter is to be used with guide wires with maximum 0.014" diameter in rapid exchange configuration.

CONTENTS

One PTCA Dilatation Catheter, one instructions for use.

MEDICAL INDICATIONS

The PTCA Dilatation Catheter is indicated for:

- coronary artery stenosis
- bypass graft stenosis
- post dilatation of balloon expandable coronary stents

INTENDED PURPOSE

The PTCA Catheter is intended for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. The balloon dilatation catheter is also indicated for post deployment expansion of balloon expandable coronary stents.

TARGET USER GROUP

PTCA Dilatation Catheters should only be used by physicians who are experienced in coronary angiography and the PTCA technique.

TARGET PATIENT POPULATION

The product shall only be used in adult persons above 18 years of age.

CONTRAINDICATIONS

The PTCA Dilatation Catheter is contraindicated for:

- non-candidates for coronary bypass graft surgery
- lesions that cannot be reached or treated with the system
- unprotected left main coronary artery
- coronary artery spasm in the absence of significant stenosis

WARNINGS

- For single patient, single procedure use only. **Do not resterilize or reuse the device.** Resterilization or reuse may compromise device performance and increase the risk of inappropriate resterilization and cross contamination.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- For patients who are not acceptable candidates for coronary artery bypass graft surgery careful consideration is required, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- To reduce the potential for vessel damage the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis. Consider both balloon diameters if kissing balloon technique is performed.
- Do not use the PTCA Dilatation Catheter if its package has been opened or damaged.
- Use the PTCA Dilatation Catheter prior to the "Use-by Date" specified on the package.
- Do not use or straighten a PTCA Dilatation Catheter if the shaft is kinked as this may lead to breakage of the shaft.
- Do not expose the PTCA Dilatation Catheter to any organic solvents, e.g. alcohol as they may have an impact on the properties of the PTCA Dilatation Catheter.
- Use only appropriate balloon inflation medium. To prevent the possibility of an air embolus, never use air or any gaseous medium to inflate the balloon.
- Do not advance or retract the PTCA Dilatation Catheter if unusual resistance is felt. When excessive force needs to be applied to move the PTCA Dilatation Catheter, stop the procedure before proceeding. Determine the root cause of the resistance before moving the catheter. Advancement by force may result in damage to the vessel and/or the PTCA Dilatation Catheter. This may necessitate recovery of fragments.
- Do not advance or retract the PTCA Dilatation Catheter unless it is fully deflated.
- Balloon pressure should not exceed the rated burst pressure (RBP) indicated on the package label for each balloon. The RBP is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over-pressurization. If inflated above the rated burst pressure, there is a risk of catheter burst.
- Do not torque the PTCA Dilatation Catheter more than one (1) full turn.
- Bench testing was conducted with a standard balloon expandable stent, e.g. Biotronik Prokinetic. Consideration should be taken when the device is used with different manufacturer's stents due to differences related to the stent design, inner surface and other characteristics. All stents should be deployed in accordance with the manufacturer's instructions for use.
- The balloon may be inflated several times, but it should not be inflated more than three times to the specified rated burst pressure. If inflated more than three times to the rated burst pressure, there is a risk of catheter burst.

PRECAUTIONS

- PTCA Dilatation Catheters should only be used by physicians who are experienced in coronary angiography and the PTCA technique.
- During the procedure, appropriate anticoagulation and antiplatelet administration must be provided to the patient. Coronary vasodilator therapy is usually provided to the patient. Antiplatelet aggregation treatment and vasodilators are usually continued in all cases for a period of time determined by the physician.
- Prior to angioplasty, all equipment to be used for the procedure, including the PTCA Dilatation Catheter, should be carefully examined to verify proper performance.
- If resistance during insertion of the PTCA Dilatation Catheter is encountered do not use force to advance the catheter as this may result in damage to the catheter.
- When the PTCA Dilatation Catheters is in the body, it should only be manipulated while under fluoroscopic observation with radiographic equipment that provides high quality images.
- If resistance during removal of the PTCA Dilatation Catheter is encountered do not use force to retract the catheter as this may result in damage to the catheter. Simultaneous removal of the entire system (e.g. guide wire and guiding catheter) might be advised.
- The safety and effectiveness of this PTCA Dilatation Catheter for the treatment of in-stent restenosis have not been established.

ADVERSE EFFECTS AND RISKS OF USE

These can include and are not be limited to:

- Acute myocardial infarction,
- Arteriovenous fistula,
- Arrhythmias including ventricular fibrillation,
- Coronary artery spasm,
- Coronary vessel dissection, perforation, rupture or other injury,
- Death,
- Drug reaction or allergic reaction to contrast medium, antiplatelets, anticoagulants,
- Embolism,
- Hemorrhage or hematoma,
- Hypo- / hypertension,
- Infection,
- Restenosis of the dilated artery,
- Stroke, air embolism and embolism due to fragmentation of thrombotic or atherosclerotic material,
- Total occlusion of the coronary artery or bypass graft,
- Unstable angina.

DIRECTIONS FOR USE

PREPARATION OF THE PTCA DILATATION CATHETER

To properly prepare the catheter for use, the following steps should be completed:

1. Remove PTCA Dilatation Catheter gently out of packaging (sterile packaging and dispenser).
2. Fill a plastic syringe with 6 to 8 ml of inflation medium. Attach the syringe to the luer lock connector located on the proximal end of the PTCA Dilatation Catheter.
3. Remove the PROTECTIVE SHEATH from the balloon. Assure that PACKAGING WIRE is also removed.
NOTE: Do not use excessive force to remove the protective sheath or packaging wire.
4. Pull a vacuum on the syringe to remove any air in the catheter lumen or in the balloon. Gentle tapping of the PTCA Dilatation Catheter may help to remove air from the balloon and the PTCA Dilatation Catheter into the syringe.
5. Flush the guide wire lumen of the PTCA Dilatation Catheter.

INSERTION TECHNIQUE

6. Place the guiding catheter, with a hemostasis valve attached, in the orifice of the target coronary artery.
7. Advance a guide wire through the guiding catheter to reach and cross the target lesion. Advance the distal tip of the PTCA Dilatation Catheter over the proximal end of the guide wire. Ensure that the guide wire exits the PTCA Dilatation Catheter through the guide wire exit port.
8. The hemostasis valve of the guiding catheter should be gradually tightened to control back flow.
NOTE: Excessive valve tightening may affect balloon inflation / deflation time as well as movement of the guide wire.
9. Track the PTCA Dilatation Catheter over the guide wire to cross the lesion using the radiopaque marker(s) to locate the balloon across the lesion under fluoroscopy.
PRECAUTION: If resistance during insertion is encountered do not use force to advance the PTCA Dilatation Catheter as this may result in damage to the PTCA Dilatation Catheter. Stop the procedure and determine the cause.
NOTE: Do not torque the PTCA Dilatation Catheter more than one (1) full turn.

BALLOON INFLATION

10. Inflate the balloon first to a low pressure to confirm location. Once proper balloon position is set under fluoroscopy, high pressure inflation may be initiated.
NOTE: Do not exceed the rated burst pressure indicated on the label.
11. After the first inflation and subsequent deflation, assess distal coronary blood flow by fluoroscopy (contrast medium injection through the guiding catheter while the balloon is deflated).
NOTE: A significant stenosis may require successive inflations. If necessary, the PTCA Dilatation Catheter may be exchanged without removing the guide wire and without an exchange wire.

CATHETER REMOVAL

12. After dilation, fully deflate the balloon and confirm that the balloon is fully deflated. Apply negative pressure until the balloon is fully deflated as observed under fluoroscopy.
13. Loosen the hemostatic valve on the guiding catheter. Hold the guide wire position and withdraw the PTCA Dilatation Catheter into the guiding catheter. Remove the PTCA Dilatation Catheter from the system.
PRECAUTION: If resistance is encountered during removal do not use force to retract the device as this may result in damage to the PTCA Dilatation Catheter. Simultaneous removal of the entire system (e.g. PTCA Dilatation Catheter, guide wire and guiding catheter) might be advised.
14. Perform a post angioplasty angiogram by contrast medium injection through the guiding catheter.
NOTE: Inspect the balloon and the PTCA Dilatation Catheter for any potential damage after wiping clean with gauze soaked with sterile normal saline.
15. Dispose of PTCA Dilatation Catheter as per internal hospital or national requirements and regulations. In case of any damage to the PTCA Dilatation Catheter of unknown cause, return the device to the manufacturer.

REFERENCES

The physician should consult recent literature on current medical practice on PTCA, such as published by the American College of Cardiology, the American Heart Association or European Society of Cardiology. Summary of Safety and Clinical Performance will be available on Eudamed.

DISCLAIMER

ALTHOUGH THE PTCA DILATATION CATHETER HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, SIS MEDICAL AG AND ITS AFFILIATES HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. SIS MEDICAL AG AND ITS AFFILIATES, THEREFORE, DISCLAIM ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PTCA DILATATION CATHETER, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. SIS MEDICAL AG AND ITS AFFILIATES SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PTCA DILATATION CATHETER, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND SIS MEDICAL AG AND ITS AFFILIATES TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PTCA DILATATION CATHETER.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Any serious incident involving this product shall be reported to the manufacturer and to the competent authority of the country where the product was used.

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SIS Medical AG
SIS Medical GmbH